

SureSmile

Great smiles in less time.

Should I be more careful with these wires?

Because your brackets work together with the wire to move your teeth precisely, please closely follow the instructions you received at the beginning of treatment to avoid accidentally damaging or knocking off your brackets. If you lose brackets after receiving your wires, you may need to return for another scan followed by a new set of wires. Not only is that inconvenient for you, but it lengthens your treatment. For additional information on braces care go to www.montanosmile.com/what-to-expect.php



Are SureSmile wires more comfortable?

Many patients have found that SureSmile wires are more comfortable for two reasons. First, the discomfort of being in braces is largely caused by wire changes. Since SureSmile wires are custom-designed, there is no need to constantly change out wires to advance treatment.

Second, the SureSmile process using a wire material called Copper Ni-Ti®. It has gained popularity for being flexible yet very effective in delivering sufficient forces to move teeth consistently. It is also heat-activated, which means it responds to body heat. Here's a tip — drink something hot if you are in a hurry to finish treatment, or drink something cold to alleviate pressure when the wires are first installed.



While Copper Ni-Ti has many advantages to recommend it for orthodontic treatment, it is almost impossible to bend reliably by hand. Due to this limitation, Copper Ni-Ti could not be used once a patient progressed past the initial stages of treatment.

With SureSmile, a robot heats the wire to high temperatures making it malleable enough to place the custom bends. Now this more comfortable wire material can be used at any point in treatment.

How do SureSmile wires shorten treatment time?

Instead of parceling out your treatment over two years of office visits, SureSmile allows Dr. Montano to plan your treatment to the end on the computer and consolidate months of appointments.

Your SureSmile custom wires achieve "a great smile in less time" because they are designed from Dr. Montano's computerized treatment plan, then bent precisely using advanced robotics.



9330 Stockdale HWY, Suite 200
Bakersfield, Ca. 93311

665-7600 or montanosmile@arrival.net
WWW.MONTANOSMILE.COM